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

Xcela Ports	Smart Port	Xcela Plus Ports	BioFlo Ports	SmartPort+ and SmartPort Plastic
				
- PASV - Vortex Technology - BioFlo Catheter - Strategic Fit	- PASV - Plastic Port Body - BioFlo Catheter - Vortex Technology	- BioFlo Catheter - Vortex Technology - Plastic Port Body - PASV	- Vortex Technology - Cost Effective - BioFlo Catheter - PASV	- PASV - BioFlo Catheter - Vortex Technology - Plastic Port Body - Cost Effective
✗ ✗ ✗ ✗	✗ ✗ ✗ ✓	✗ ✗ ✓ ✓	✗ ✗ ✓ ✓	✗ ✓ ✓ ✓ ✓

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




SmartPort PLASTIC

Port Body	Port Design	Catheter Size	Catheter Material
Plastic	-	6F, 8F	Polyurethane

All configurations come available in the following options:
 Catheter: Attached or Detached
 Suture holes: Filled or Non-Filled

Catheter Size	Introducer Size
5FR	6FR
6FR	6FR
8FR	8FR

Kit


SmartPort+ and Plastic Kit

- Port Body
- Catheter
- 18 Gauge Introducer needle
- 0.038" J-tipped guidewire
- Greatbatch valved sheath
- 9" Metal Tunneler
- Plastic cannula (blunt needle)
- 12 mL luer lock syringe
- 12 mL slip lock syringe
- 22 gauge straight huber needle
- 22 gauge 90deg huber needle
- Vein pick

What's New?

- 5FR Catheter
- Locking Collars

Estudios

Toledo Study

Title: A Randomized, Prospective Trial of Conventional Vascular Ports Vs. The Vortex “Clear-Flow” Reservoir Port in Adult Oncology Patients

Authors: Barbara Stevens, RN; Sue Ellen Barton, RN; Marjorie Brechbill, RN; Sue Moenter, RN; Anna Lou Piel, RN; Darcel Shankle, RN

Where: Toledo Oncology Clinic, Toledo, OH

What: 95 Patients: **48 received Vortex port**, 47 received a conventional port. 18,144 dwell days studied.

Findings: Occlusions were reduced by 73 %. The need for interventions was reduced by 69 %.

The following article appeared in the *Journal of Vascular Access Devices (JVAD)*, the official publication of the National Association of Vascular Access Networks (NAVAN). In September 2003, the name of *JVAD* was changed to the *Journal of the Association for Vascular Access* and the name of NAVAN was changed to the Association for Vascular Access (AVA). For more information, contact AVA at 1-877-921-3821 or visit our website at www.avanet.org.

A Randomized, Prospective Trial of Conventional Vascular Ports vs. The Vortex “Clear-Flow” Reservoir Port in Adult Oncology Patients


Barbara Stevens, RN, OCN, CCRA
Sue Ellen Barton, RN, OCN
Marjorie Brechbill, RN, OCN

Sue Moenter, RN, OCN
Anna Lou Piel, RN
Darcel Shankle, RN, OCN

Introduction
Vascular access ports were developed to provide long-term intravenous access with limited peripheral access, combined with the need for frequent venipuncture, in oncology patients receiving long-term, intensive therapy. Since their introduction in

1983, implantable vascular access ports have become a standard of use in the treatment of long-term venous therapy. Over the past decade, the use of implantable ports has increased significantly. The literature does reveal many problems associated with their use. The most common complications associated with frequent venipuncture, in oncology patients receiving long-term, intensive therapy, are occlusions and thrombus formations. Vascular access ports have proven to be very beneficial in overcoming the


need for frequent venipuncture in the treatment of long-term venous therapy. Over the past decade, the use of implantable ports has increased significantly. The literature does reveal many problems associated with their use. The most common complications associated with frequent venipuncture, in oncology patients receiving long-term, intensive therapy, are occlusions and thrombus formations. Vascular access ports have proven to be very beneficial in overcoming the

Call Points

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
Quirófano


Radiología

Oncología

Investigating Clinical Call Points

Competencia

Bard Ports

Products

- Power Port
- Power Port Slim
- Power Port Duo
- Power Port VUE
- Vas-Cath Port (value model)
- Vaccess (Low cost plastic)
- Clearvue (slim)

Strengths

- Palpation Nubs
- Only Silicone Port body on market
- Full offering of all products
- GPO and IDN Contracts
- Large Sales and Corporate Accounts Team
- Multi-Billion Dollar Marketing Machine
- 90 Day flush

Weaknesses

- No Anti-Thrombogenic Technology
- No Proximal Valve Port Technology
- Palpation Nub
- Medical Grade Plastic

Braun

Products

- Celsite PSU – plastic
- Celsite EPOXY – titanium
- Celsite Discreet – titanium
- Celsite Valved** valved port

Strengths

- Name recognition
- Breadth of offering
- Valved options
- ECG tip location offering
- Arterial offering

Weaknesses

- No Anti-Thrombogenic Technology
- No Vortex Technology

MedComp Ports

Products

- ProFuse CT
- Dignity

Strengths

- Price
- Low Profile
- Dignity is a solid plastic port
- Small Sales Bag
- Strong in NE region
- Pediatric indication

Weaknesses

- No Anti-Thrombogenic Ports
- Not on many GPO/IDN Contracts
- No Valved Port options
- Regionally strong, not nationally
- Large Territories; smaller sales force
- Sales Force Turnover High

Smith's Medical Ports

Products

- PORT-A-CATH
- PORT-A-CATH II
- P.A.S. Port

Strengths

- Peripheral Port Offering
- Ultra-Lock Locking Collar
- Hybrid Port Offering
- Focused in OR
- Small Sales Bag

Weaknesses

- No Anti-Thrombogenic Ports
- Not on many GPO/IDN Contracts
- No Valved Port options
- Regionally strong, not nationally
- Large Territories; Smaller sales force