

YOUR CLASSIC CVC RANGE



REDUCE
PATIENT **RISKS**
WITH ALTIUS



A CLASSIC CATHETER WITH A SUPERIOR INSERTION KIT

Altius Classic range provides all the benefits of a traditional catheter with a kit designed to reduce risks during the catheter insertion.

Altius Classic is manufactured using our proprietary Multitube manufacturing technique.

The Multitube technique allows to manufacture perfectly round lumens and therefore achieving laminar flow.

This allows for higher flows than conventional extruded catheters specially in 4 & 5 lumen configuration.

HIGH BARIUM CONTENT

Altius Classic is manufactured with a high barium content to aid in the visualisation of the catheter during X ray tip confirmation.



MULTITUBE TECHNOLOGY

upto 21% higher flows than traditional 5 lumen extruded catheters.*



THE RANGE

1 to 5 lumens with different lengths.
See the Ordering Guide for Offering & Flow rates.

*comparison of market leader 5 lumen 8.5 french 16cm published data with Altius Classic 5 lumen 8.5Fr.

SUPERIOR INSERTION KIT

Altius Classic superior insertion kit has been designed to reduce the risks during insertion.



Super Sharp Echogenic Needle: Easy visualisation under the ultrasound to prevent posterior wall puncture.



5m Guidewire Introducer: Luerslip syringe for closed system.



Altius Guidewire Dispenser: Ergonomic design to facilitate guidewire advancement even with wet gloves.



Kink Resistant Nitinol Guidewire: To prevent guidewire kinking even in difficult patient's anatomy.



N11 Mini Scalpel: To nick skin prior to dilation.



Hydrophilic Dilator: For smooth dilation to prevent insertion trauma.



Secondary Fixation Wing

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

MIDDLE EAST, ASIA & AFRICA

Tel: +97 16 552 8146

Fax: +97 16 552 8147

Email: mea@kimal.co.uk

INTERNATIONAL

Tel: +44 (0)845 437 9543

Fax: +44 (0)845 437 9544

Email: internationalenquiries@kimal.co.uk

© 2021 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC.

Kimal PLC reserve the right to change the specification of its products at any time.
Kimal PLC, Arundel Road, Uxbridge, UB8 2SA.

??????? | NOVEMBER 2021

KIMAL
delivering healthcare innovation