



Neopic

Neonatal Peripherally Inserted Central Catheter
(Neonatal PICC)

GAMAMED

NEOPIC

Neonatal Peripherally Inserted Central Catheter (Neonatal PICC)

- **Prevents Repeated Punctures with IV Cannulas:** Preserves the integrity of peripheral veins.
- **Suitable for the Infusion of:** Parenteral nutrition, hyperosmolar solutions, and vasoactive drugs.
- **Hemocompatible:** Ensures safe interaction with blood components.
- **Soft, Well-Tolerated Catheter:** Helps prevent phlebitis and enhances patient comfort.
- **Tensile Strength Exceeds ISO 10555-3 Standards:** Reduces the risk of catheter breakage.



CONTENTS

- 1 Catheter.
- 1 Guidewire.
- 1 Dilator.
- 1 Introducer Needle.
- 1 Blue Introducer syringe.
- 1 Injection Cap.
- 1 Syringe
- 1 Injection Needle.

NEOPIC (Neonatal PICC)

CODE	DESCRIPTION
NEOP-19140	PICC Cateter (Pur) 1.9Fr 1LUMEN 40cm
NEOP-20140	PICC Cateter (Pur) 2Fr 1 LUMEN 40cm
NEOP-20240	PICC Cateter (Pur) 2Fr 2 LUMEN 40cm

Gamamed Vascular Systems Inc

Redefining Vascular Wellness

Operational Headquarters

815 14th Street SW, Suite C250,
Loveland CO 80537

Email: info@gamavascular.com

phone: +1 (307) 998 - 9732

www.gamavascular.com

From intravenous catheters to central venous access systems, our solutions are trusted by hospitals, clinics, and healthcare providers around the world.