




Precision Flow Long Term

Hemodialysis Catheter Kit





Gamamed Vascular Systems Inc
Redefining Vascular Wellness

PRECISION FLOW LONG TERM HEMODIALYSIS CATHETER KIT

SYMMETRICAL TIP

- > **Symmetrical Tip Design:** Promotes catheter patency and optimal flow dynamics.
- > **High-Strength Material with Optimized Inner Lumen Design:** Combines a durable, high-tensile-strength material with an engineered lumen shape to maximize flow rate without compromising flexibility or kink resistance.
- > **Minimal Recirculation:** Ensures efficient dialysis performance.
- > **PU Tubing:** Durable and biocompatible polyurethane construction.



Symmetrical-tip

Hemodialysis Hemo-Symmetrics	
CODE	DESCRIPTION
SYM-008-2-11-28	Chronic Hemodialysis Catheter 8Fr X 2lumen X 11cm X 28cm long. Tip to Cetf
SYM-008-2-13-30	Chronic Hemodialysis Catheter 8Fr X 2lumen X 13cm X 30cm long. Tip to Cetf
SYM-008-2-15-32	Chronic Hemodialysis Catheter 8Fr X 2lumen X 15cm X 32cm long. Tip to Cetf
SYM-010-2-12-29	Chronic Hemodialysis Catheter 10Fr X 2lumen X 12cm X 29cm long. Tip to Cetf
SYM-010-2-13-30	Chronic Hemodialysis Catheter 10Fr X 2lumen X 13cm X 30cm long. Tip to Cetf
SYM-010-2-15-32	Chronic Hemodialysis Catheter 10Fr X 2lumen X 15cm X 32cm long. Tip to Cetf
SYM-010-2-16-33	Chronic Hemodialysis 10Fr X 2lumen X 16cm X 33cm long. Tip to Cetf
SYM-012-2-16-33	Chronic Hemodialysis 12Fr X 2lumen X 16cm X 33cm long. Tip to Cetf
SYM-012-2-19-36	Chronic Hemodialysis Catheter 12Fr X 2lumen X 19cm X 36cm long. Tip to Cetf
SYM-012-2-23-40	Chronic Hemodialysis Catheter 12Fr X 2lumen X 23cm X 40cm long. Tip to Cetf
SYM-135-2-16-33	Chronic Hemodialysis Catheter 13.5Fr X 2lumen X 16cm X 33cm long. Tip to Cetf
SYM-135-2-19-36	Chronic Hemodialysis Catheter 13.5Fr X 2lumen X 19cm X 36cm long. Punta a Cetf
SYM-135-2-23-40	Chronic Hemodialysis Catheter 13.5Fr X 2lumen X 23cm X 40cm long. Tip to Cetf
SYM-135-2-27-44	Chronic Hemodialysis Catheter 13.5Fr X 2lumen X 27cm X 44cm long. Tip to Cetf
SYM-135-2-28-45	Chronic Hemodialysis Catheter 13.5Fr X 2lumen X 28cm X 45cm long. Tip to Cetf



Gamamed Vascular Systems Inc
Redefining Vascular Wellness

Precision Flow


Hemodialysis Catheter Kit

Hemodialysis Hemo-Symmetrics	
CODE	DESCRIPTION
SYM-135-2-33-50	Chronic Hemodialysis Catheter 13.5Fr X 2lumen X 33cm X 50cm long. Tip to Cetf
SYM-145-2-19-36	Chronic Hemodialysis Catheter 14.5Fr X 2lumen X 19cm X 36cm long. Tip to Cetf
SYM-145-2-23-40	Chronic Hemodialysis Catheter 14.5Fr X 2lumen X 23cm X 40cm long. Tip to Cetf
SYM-145-2-28-45	Chronic Hemodialysis Catheter 14.5Fr X 2lumen X 28cm X 45cm long. Tip to Cetf
SYM-145-2-33-50	Chronic Hemodialysis Catheter 14.5Fr X 2lumen X 33cm X 50cm long. Tip to Cetf
SYM-145-2-44-61	Chronic Hemodialysis Catheter 14.5Fr X 2lumen X 44cm X 61cm long. Tip to Cetf
SYM-145-2-55-72	Chronic Hemodialysis Catheter 14.5Fr X 2lumen X 55cm X 72cm long. Tip to Cetf
SYM-155-2-19-36	Chronic Hemodialysis Catheter 15.5Fr X 2lumen X 19cm X 36cm long. Tip to Cetf
SYM-155-2-23-40	Chronic Hemodialysis Catheter 15.5Fr X 2lumen X 23cm X 40cm long. Tip to Cetf
SYM-155-2-28-45	Chronic Hemodialysis Catheter 15.5Fr X 2lumen X 28cm X 45cm long. Tip to Cetf
SYM-155-2-33-50	Chronic Hemodialysis Catheter 15.5Fr X 2lumen X 33cm X 50cm long. Tip to Cetf
SYM-155-2-44-61	Chronic Hemodialysis Catheter 15.5Fr X 2lumen X 44cm X 61cm long. Tip to Cetf
SYM-155-2-55-72	Chronic Hemodialysis Catheter 15.5Fr X 2lumen X 55cm X 72cm long. Tip to Cetf




STEP-TIP

- **Soft, Comfortable Silicone Tubing:**
Ensures patient comfort and long-term biocompatibility.
- **Step-Tip Design:**
Reduces the risk of blood clot formation.
- **PU Tubing:**
Durable and biocompatible polyurethane construction.



Step-tip

Hemodialysis Step-Tip	
CODE	DESCRIPTION
NEF-145-2-19-36	Chronic Hemodialysis Catheter 14.5Fr 2 Lumen 19cm X 36cm long.
NEF-145-2-23-40	Chronic Hemodialysis Catheter 14.5Fr 2 Lumen 23cm X 40cm long.
NEF-145-2-28-45	Chronic Hemodialysis Catheter 14.5Fr 2 Lumen 28cm X 45cm long.
NEF-145-2-33-50	Chronic Hemodialysis Catheter 14.5Fr 2 Lumen 33cm X 50cm long.
NEF-145-2-44-61	Chronic Hemodialysis Catheter 14.5Fr 2 Lumen 44cm X 61cm long.
NEF-145-2-55-72	Chronic Hemodialysis Catheter 14.5Fr 2 Lumen 55cm X 72cm long.





Gamamed Vascular Systems Inc
Redefining Vascular Wellness

Precision Flow

Hemodialysis Catheter Kit

Here's a rephrased version:

Our long-term hemodialysis kits, featuring both symmetric tip and step tip options, contain all the essential components required by medical staff for patient insertion

CONTENTS

1 Long-term hemodialysis catheter, symmetrical tip/Step Tip.
1 10 Ga × 7 cm introducer needle.
1 Injection needle.
1 16 Fr × 13 cm peel-off sheath with valve.
1 Tunneling device.
1 0.038" × 76 cm Nitinol guidewire with straight/J tip.

1 10cc syringe.
1 Raulerson syringe.
1 wound dressing.
2 heparin caps.
2 end caps.



Gamamed Vascular Systems Inc
Redefining Vascular Wellness



Operational Headquarters

815 14th Street SW, Suite C250,
Loveland CO 80537

Email: info@gamavascular.com

phone: +1 307 998 97 21

www.gamavascular.com

From intravenous catheters to central venous access systems, our solutions are trusted by hospitals, clinics, and healthcare providers around the world.