



Altius Classic Central Venous Catheter



Product Code	Description	Lumens	French size	Length (cm)	Priming Volume (ml)					Flow Rates (ml/hr)				
					Proximal	Distal	Medial 1	Medial 2	Medial 3	Proximal	Distal	Medial 1	Medial 2	Medial 3
KCS115702	Altius Classic 2 Lumen CVC Set 7Fr 15cm	2	7	15	0.35	0.40				1600	8500			
KCS120702	Altius Classic 2 Lumen CVC Set 7Fr 20cm	2	7	20	0.42	0.45				1350	7200			
KCS130702	Altius Classic 2 Lumen CVC Set 7Fr 30cm	2	7	30	0.50	0.53				1100	6250			
KCS115802	Altius Classic 2 Lumen CVC Set 8Fr 15cm	2	8	15	0.38	0.43				3200	8500			
KCS120802	Altius Classic 2 Lumen CVC Set 8Fr 20cm	2	8	20	0.45	0.50				2750	7200			
KCS115703	Altius Classic 3 Lumen CVC Set 7Fr 15cm	3	7	15	0.24	0.31	0.27			1600	5800	1600		
KCS120703	Altius Classic 3 Lumen CVC Set 7Fr 20cm	3	7	20	0.30	0.38	0.29			1350	5100	1350		
KCS130703	Altius Classic 3 Lumen CVC Set 7Fr 30cm	3	7	30	0.35	0.45	0.35			1100	4300	1100		
KCS115853	Altius Classic 3 Lumen CVC Set 8.5Fr 15cm	3	8.5	15	0.26	0.35	0.28			3200	10200	3200		
KCS120853	Altius Classic 3 Lumen CVC set 8.5Fr 20cm	3	8.5	20	0.32	0.40	0.33			2750	8500	2750		
KCS115854	Altius Classic 4 Lumen CVC Set 8.5Fr 15cm	4	8.5	15	0.28	0.45	0.35	0.28		1600	10200	2,500	1600	
KCS120854	Altius Classic 4 Lumen CVC Set 8.5Fr 20cm	4	8.5	20	0.32	0.53	0.40	0.32		1350	8500	1700	1350	
KCS130854	Altius Classic 4 Lumen CVC Set 8.5Fr 30cm	4	8.5	30	0.40	0.60	0.48	0.40		1100	8100	1350	1100	
KCS115855	Altius Classic 5 Lumen CVC Set 8.5Fr 15cm	5	8.5	15	0.25	0.35	0.28	0.26	0.25	780	10500	1250	1250	780
KCS120855	Altius Classic 5 Lumen CVC Set 8.5Fr 20cm	5	8.5	20	0.28	0.40	0.32	0.30	0.28	600	8500	1050	1050	600
KCS130855	Altius Classic 5 Lumen CVC Set 8.5Fr 30cm	5	8.5	30	0.33	0.48	0.40	0.36	0.33	550	8100	850	850	550

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

© 2024 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC.
Kimal PLC reserve the right to change the specification of its products at any time.
Kimal PLC, Arundel Road, Uxbridge, UB8 2SA.

UK
Tel: +44 (0)845 437 9542
Fax: +44 (0)845 437 9544
Email: customer care@kimal.com

International
Tel: +44 (0)845 437 9543
Fax: +44 (0)845 437 9549
Email: internationalenquiries@kimal.com