



Altius ProActiv+ High Pressure Central Venous Catheter



Product Code	Description	Lumens	French size	Length (cm)	Priming Volume (ml)						Flow Rates (ml/hr)					
					Proximal	Distal	Medial 1	Medial 2	Medial 3	Medial 4	Proximal	Distal	Medial 1	Medial 2	Medial 3	Medial 4
K2CS215501P	ALTIUS ProActiv+ HP 1 Lumen CVC Set 5Fr 15cm	1	5	15		0.3						3000(16G)				
K2CS220501P	ALTIUS ProActiv+ HP 1 Lumen CVC Set 5Fr 20cm	1	5	20		0.3						2800(16G)				
K2CS215701P	ALTIUS ProActiv+ HP 1 Lumen CVC Set 7Fr 15cm	1	7	15		0.55						9300(14G)				
K2CS220701P	ALTIUS ProActiv+ HP 1 Lumen CVC Set 7Fr 20cm	1	7	20		0.6						7500(14G)				
K2CS215702P	ALTIUS ProActiv+ HP 2 Lumen CVC Set 7Fr 15cm	2	7	15	0.3	0.35					2500(17G)	5800(14G)				
K2CS220702P	ALTIUS ProActiv+ HP 2 Lumen CVC Set 7Fr 20cm	2	7	20	0.3	0.35					2050(17G)	5100(14G)				
K2CS215802P	ALTIUS ProActiv+ HP 2 Lumen CVC Set 8Fr 15cm	2	8	15	0.3	0.4					3200(16G)	6600(14G)				
K2CS220802P	ALTIUS ProActiv+ HP 2 Lumen CVC Set 8Fr 20cm	2	8	20	0.3	0.4					2700(16G)	6100(14G)				
K2CS215703P	ALTIUS ProActiv+ HP 3 Lumen CVC Set 7Fr 15cm	3	7	15	0.25	0.3	0.25				1600(18G)	3200(16G)	1600(18G)			
K2CS215803P	ALTIUS ProActiv+ HP 3 Lumen CVC Set 8Fr 15cm	3	8	15	0.25	0.25	0.3				1600(18G)	6600(14G)	1600(18G)			
K2CS220703P	ALTIUS ProActiv+ HP 3 Lumen CVC Set 7Fr 20cm	3	7	20	0.25	0.3	0.25				1350(18G)	2800(16G)	1350(18G)			
K2CS220803P	ALTIUS ProActiv+ HP 3 Lumen CVC Set 8Fr 20cm	3	8	20	0.25	0.25	0.3				1350(18G)	5100(14G)	1350(18G)			
K2CS230703P	ALTIUS ProActiv+ HP 3 Lumen CVC Set 7Fr 30cm	3	7	30	0.25	0.3	0.25				1100(18G)	2350(16G)	1100(18G)			
K2CS215904P	ALTIUS ProActiv+ HP 4 Lumen CVC Set 9Fr 15cm	4	9	15	0.25	0.25	0.25	0.25			1600(17G)	6600(14G)	1600(17G)	1600(17G)		
K2CS220904P	ALTIUS ProActiv+ HP 4 Lumen CVC Set 9Fr 20cm	4	9	20	0.25	0.25	0.25	0.25			1350(17G)	6050(14G)	1350(17G)	1350(17G)		
K2CS215905P	ALTIUS ProActiv+ HP 5 Lumen CVC Set 9Fr 15cm	5	9	15	0.2	0.3	0.2	0.2	0.2		1250(18G)	3200(16G)	1250(18G)	1250(18G)	1250(18G)	
K2CS220905P	ALTIUS ProActiv+ HP 5 Lumen CVC Set 9Fr 20cm	5	9	20	0.2	0.3	0.2	0.2	0.2		1100(18G)	2800(16G)	1100(18G)	1100(18G)	1100(18G)	
K2CS230905P	ALTIUS ProActiv+ HP 5 Lumen CVC Set 9Fr 30cm	5	9	30	0.2	0.3	0.2	0.2	0.2		840(18G)	2200(16G)	840(18G)	840(18G)	840(18G)	

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

© 2024 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC.
Kimal PLC reserve the right to change the specification of its products at any time.
Kimal PLC, Arundel Road, Uxbridge, UB8 2SA.

UK
Tel: +44 (0)845 437 9542
Fax: +44 (0)845 437 9544
Email: customer care@kimal.com

International
Tel: +44 (0)845 437 9543
Fax: +44 (0)845 437 9549
Email: internationalenquiries@kimal.com