

OVER DILATOR

Product Code	Description	Lumens	French size	Length (cm)	Priming Volume overall catheter (ml)			Flow Rates (ml/hr)		
					Proximal	Medial	Distal	Proximal	Medial	Distal
KP-CTSL-0355-KOD	Kimal PICC CT Compatible 3Fr x 55cm Single Lumen Kit Over Dilator	1	3	55			0.40			400
KP-CTSL-0455-KOD	Kimal PICC CT Compatible 4Fr x 55cm Single Lumen Kit Over Dilator	1	4	55			0.90			800
KP-CTSL-0555-KOD	Kimal PICC CT Compatible 5Fr x 55cm Single Lumen Kit Over Dilator	1	5	55			0.90			1220
KP-CTDL-0555-KOD	Kimal PICC CT Compatible 5Fr x 55cm Dual Lumen Kit Over Dilator	2	5	55	0.57		0.57	220		680
KP-CTDL-0655-KOD	Kimal PICC CT Compatible 6Fr x 55cm Dual Lumen Kit Over Dilator	2	6	55	0.73		0.73	450		1100
KP-CTTL-0655-KOD	Kimal PICC CT Compatible 6Fr x 55cm Triple Lumen Kit Over Dilator	3	6	55	0.20	0.20	0.65	170	170	900

OVER NEEDLE

KP-CTSL-0355-KON	Kimal PICC CT Compatible 3Fr x 55cm Single Lumen Kit Over Needle	1	3	55			0.40			400
KP-CTSL-0455-KON	Kimal PICC CT Compatible 4Fr x 55cm Single Lumen Kit Over Needle	1	4	55			0.90			800
KP-CTSL-0555-KON	Kimal PICC CT Compatible 5Fr x 55cm Single Lumen Kit Over Needle	1	5	55			0.90			1220
KP-CTDL-0555-KON	Kimal PICC CT Compatible 5Fr x 55cm Dual Lumen Kit Over Needle	2	5	55	0.63		0.63	220		680
KP-CTDL-0655-KON	Kimal PICC CT Compatible 6Fr x 55cm Dual Lumen Kit Over Needle	2	6	55	0.73		0.73	450		1100
KP-CTTL-0655-KON	Kimal PICC CT Compatible 6Fr x 55cm Triple Lumen Kit Over Needle	3	6	55	0.20	0.20	0.65	170	170	900

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

© 2024 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC.
Kimal PLC reserve the right to change the specification of its products at any time.
Kimal PLC, Arundel Road, Uxbridge, UB8 2SA.

UK
Tel: +44 (0)845 437 9542
Fax: +44 (0)845 437 9544
Email: customer-care@kimal.com

International
Tel: +44 (0)845 437 9543
Fax: +44 (0)845 437 9549
Email: internationalenquiries@kimal.com