

ADULT / SINGLE LUMEN / STRAIGHT DESIGN				FLOW RATES										Priming Volume		
				100ml / min		150 ml / min		200ml / min		300ml / min		400ml / min				
Kimal Product Code	Description	French	Length (cm)	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Medial	Venous
KFE-TSL-0815-K	Short Term Single Lumen Straight Haemodialysis Catheter Kit Size 8Fr Length 15cm	8	15	-50		-70		-100						0.70		
KFE-TSL-0820-K	Short Term Single Lumen Straight Haemodialysis Catheter Kit Size 8Fr Length 20cm	8	20	-50		-70		-100						0.80		

ADULT / DUAL LUMEN / STRAIGHT DESIGN				FLOW RATES										Priming Volume		
				100ml / min		150 ml / min		200ml / min		300ml / min		400ml / min				
Kimal Product Code	Description	French	Length (cm)	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Medial	Venous
KFE-TDL-1015-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 10Fr Length 15cm	10	15					-100	120	-190	210	-185	290	0.90		1.00
KFE-TDL-1212-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 12Fr Length 12cm	12	12					-60	70	-145	160	-180	190	0.90		1.00
KFE-TDL-1215-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 12Fr Length 15cm	12	15					-50	80	-130	135	-180	190	1.00		1.10
KFE-TDL-1217-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 12Fr Length 17cm	12	17					-70	70	-165	180	-200	210	1.10		1.10
KFE-TDL-1220-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 12Fr Length 20cm	12	20					-70	70	-175	190	-210	220	1.10		1.20
KFE-TDL-1224-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 12Fr Length 24cm	12	24					-95	100	-160	155	-235	260	1.20		1.30
KFE-TDL-1415-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 14Fr Length 15cm	14	15					-50	60	-130	135	-170	180	1.10		1.10
KFE-TDL-1417-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 14Fr Length 17cm	14	17					-70	70	-145	150	-160	170	1.10		1.20
KFE-TDL-1420-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 14Fr Length 20cm	14	20					-60	60	-150	160	-200	210	1.20		1.30
KFE-TDL-1424-K	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 14Fr Length 24cm	14	24					-60	60	-150	160	-200	210	1.30		1.40

ADULT / DUAL LUMEN / CURVED DESIGN				FLOW RATES										Priming Volume		
				100ml / min		150 ml / min		200ml / min		300ml / min		400ml / min				
Kimal Product Code	Description	French	Length (cm)	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Medial	Venous
KFE-TDL-1015-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 10Fr Length 15cm	10	15					-100	120	-190	210	-285	290	0.90		1.00
KFE-TDL-1212-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 12Fr Length 12cm	12	12					-60	70	-145	160	-180	190	0.90		1.00
KFE-TDL-1215-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 12Fr Length 15cm	12	15					-50	80	-130	135	-180	190	1.00		1.10
KFE-TDL-1217-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 12Fr Length 17cm	12	17					-70	70	-165	180	-200	210	1.10		1.10
KFE-TDL-1220-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 12Fr Length 20cm	12	20					-70	70	-175	190	-210	220	1.10		1.20
KFE-TDL-1224-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 12Fr Length 24cm	12	24					-95	100	-160	155	-235	260	1.20		1.30
KFE-TDL-1415-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 14Fr Length 15cm	14	15					-50	60	-130	135	-170	180	1.10		1.10
KFE-TDL-1417-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 14Fr Length 17cm	14	17					-70	70	-145	150	-160	170	1.10		1.20
KFE-TDL-1420-JK	Short Term Dual Lumen Curved Ext Haemodialysis Catheter Kit Size 14Fr Length 20cm	14	20					-60	60	-150	160	-200	210	1.20		1.30

ADULT / DUAL LUMEN / PRE CURVED DESIGN				FLOW RATES										Priming Volume		
				100ml / min		150 ml / min		200ml / min		300ml / min		400ml / min				
Kimal Product Code	Description	French	Length (cm)	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Medial	Venous
KFE-TDL-1212-PCK	Short Term Dual Lumen Precurved Haemodialysis Catheter Kit Size 12Fr Length 12cm	12	12					-70	80	-130	130	-200	210	1.00		1.10
KFE-TDL-1215-PCK	Short Term Dual Lumen Precurved Haemodialysis Catheter Kit Size 12Fr Length 15cm	12	15					-90	95	-145	140	-210	220	1.10		1.20
KFE-TDL-1217-PCK	Short Term Dual Lumen Precurved Haemodialysis Catheter Kit Size 12Fr Length 17cm	12	17					-90	100	-150	145	-230	240	1.20		1.20
KFE-TDL-1220-PCK	Short Term Dual Lumen Precurved Haemodialysis Catheter Kit Size 12Fr Length 20cm	12	20					-95	100	-160	155	-235	260	1.20		1.30
KFE-TDL-1224-PCK	Short Term Dual Lumen Pre Curved Haemodialysis Catheter Kit Size 12Fr Length 24cm	12	24					-100	110	-170	166	-245	280	1.30		1.40
KFE-TDL-1415-PCK	Short Term Dual Lumen Precurved Haemodialysis Catheter Kit Size 14Fr Length 15cm	14	15					-60	60	-135	135	-200	210	1.20		1.20
KFE-TDL-1417-PCK	Short Term Dual Lumen Precurved Haemodialysis Catheter Kit Size 14Fr Length 17cm	14	17					-70	70	-140	140	-210	220	1.20		1.30
KFE-TDL-1420-PCK	Short Term Dual Lumen Precurved Haemodialysis Catheter Kit Size 14Fr Length 20cm	14	20					-70	70	-140	145	-220	230	1.30		1.40
KFE-TDL-1424-PCK	Short Term Dual Lumen Precurved Haemodialysis Catheter Kit Size 14Fr Length 24cm	14	24					-70	70	-140	145	-220	230	1.30		1.40

ADULT / TRIPLE LUMEN / STRAIGHT DESIGN				FLOW RATES										Priming Volume		
				100ml / min		150 ml / min		200ml / min		300ml / min		400ml / min				
Kimal Product Code	Description	French	Length (cm)	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Medial	Venous
KFE-TTL-1315-K	Short Term Triple Lumen Straight Haemodialysis Catheter Kit Size 13Fr Length 15cm	13	15					-90	90	-110	125	-160	170	1.10	0.50	1.00
KFE-TTL-1317-K	Short Term Triple Lumen Straight Haemodialysis Catheter Kit Size 13Fr Length 17cm	13	17					-90	90	-125	130	-160	170	1.15	0.55	1.05
KFE-TTL-1320-K	Short Term Triple Lumen Straight Haemodialysis Catheter Kit Size 13Fr Length 20cm	13	20					-90	90	-120	135	-170	180	1.20	0.60	1.10
KFE-TTL-1324-K	Short Term Triple Lumen Straight Haemodialysis Catheter Kit Size 13Fr Length 24cm	13	24					-100	100	-130	145	-160	160	1.25	0.65	1.15
KFE-TTL-1415-K	Short Term Triple Lumen Straight Haemodialysis Catheter Kit Size 14Fr Length 15cm	14	15					-90	90	-110	125	-160	170	1.10	0.50	1.00
KFE-TTL-1417-K	Short Term Triple Lumen Straight Haemodialysis Catheter Kit Size 14Fr Length 17cm	14	17					-90	90	-125	130	-160	170	1.15	0.55	1.05
KFE-TTL-1420-K	Short Term Triple Lumen Straight Haemodialysis Catheter Kit Size 14Fr Length 20cm	14	20					-90	90	-120	135	-170	180	1.20	0.60	1.10
KFE-TTL-1424-K	Short Term Triple Lumen Straight Haemodialysis Catheter Kit Size 14Fr Length 24cm	14	24					-100	100	-130	145	-160	160	1.25	0.65	1.15

ADULT / TRIPLE LUMEN / CURVED DESIGN				FLOW RATES										Priming Volume		
				100ml / min		150 ml / min		200ml / min		300ml / min		400ml / min				
Kimal Product Code	Description	French	Length (cm)	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Medial	Venous
KFE-TTL-1315-JK	Short Term Triple Lumen Curved Ext Haemodialysis Catheter Kit Size 13Fr Length 15cm	13	15					-90	90	-110	125	-160	170	1.10	0.50	1.10
KFE-TTL-1317-JK	Short Term Triple Lumen Curved Ext Haemodialysis Catheter Kit Size 13Fr Length 17cm	13	17					-90	90	-125	130	-160	170	1.15	0.55	1.10
KFE-TTL-1320-JK	Short Term Triple Lumen Curved Ext Haemodialysis Catheter Kit Size 13Fr Length 20cm	13	20					-90	90	-120	135	-170	180	1.20	0.60	1.15
KFE-TTL-1324-JK	Short Term Triple Lumen Curved Ext Haemodialysis Catheter Kit Size 13Fr Length 24cm	13	24					-100	100	-130	145	-160	160	1.25	0.65	1.20
KFE-TTL-1415-JK	Short Term Triple Lumen Curved Ext Haemodialysis Catheter Kit Size 14Fr Length 15cm	14	15					-90	90	-110	125	-160	170	1.10	0.50	1.00
KFE-TTL-1417-JK	Short Term Triple Lumen Curved Ext Haemodialysis Catheter Kit Size 14Fr Length 17cm	14	17					-90	90	-125	130	-160	170	1.15	0.55	1.05
KFE-TTL-1420-JK	Short Term Triple Lumen Curved Ext Haemodialysis Catheter Kit Size 14Fr Length 20cm	14	20					-90	90	-120	135	-170	180	1.20	0.60	1.10
KFE-TTL-1424-JK	Short Term Triple Lumen Curved Ext Haemodialysis Catheter Kit Size 14Fr Length 24cm	14	24					-100	100	-130	145	-160	160	1.25	0.65	1.15

ADULT / TRIPLE LUMEN / PRE CURVED DESIGN				FLOW RATES										Priming Volume		
				100ml / min		150 ml / min		200ml / min		300ml / min		400ml / min				
Kimal Product Code	Description	French	Length (cm)	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Medial	Venous
KFE-TTL-1315-PCK	Short Term Triple Lumen Pre Curved Haemodialysis Catheter Kit Size 13Fr Length 15cm	13	15					-90	90	-120	135	-170	180	1.20	0.60	1.00
KFE-TTL-1317-PCK	Short Term Triple Lumen Pre Curved Haemodialysis Catheter Kit Size 13Fr Length 17cm	13	17					-100	100	-130	145	-150	150	1.20	0.60	1.05
KFE-TTL-1320-PCK	Short Term Triple Lumen Pre Curved Haemodialysis Catheter Kit Size 13Fr Length 20cm	13	20					-100	100	-130	145	-160	160	1.25	0.65	1.10
KFE-TTL-1324-PCK	Short Term Triple Lumen Pre Curved Haemodialysis Catheter Kit Size 13Fr Length 24cm	13	24					-100	100	-130	145	-160	160	1.30	0.70	1.15
KFE-TTL-1415-PCK	Short Term Triple Lumen Pre Curved Haemodialysis Catheter Kit Size 14Fr Length 15cm	14	15					-90	90	-120	135	-170	180	1.20	0.60	1.10
KFE-TTL-1417-PCK	Short Term Triple Lumen Pre Curved Haemodialysis Catheter Kit Size 14Fr Length 17cm	14	17					-100	100	-130	145	-150	150	1.20	0.60	1.10
KFE-TTL-1420-PCK	Short Term Triple Lumen Pre Curved Haemodialysis Catheter Kit Size 14Fr Length 20cm	14	20					-100	100	-130	145	-160	160	1.25	0.65	1.15
KFE-TTL-1424-PCK	Short Term Triple Lumen Pre Curved Haemodialysis Catheter Kit Size 14Fr Length 24cm	14	24					-100	100	-130	145	-160	160	1.30	0.70	1.20

PAEDIATRIC / DUAL LUMEN / STRAIGHT DESIGN						FLOW RATES						Priming Volume		
						100ml / min		150 ml / min		200ml / min				
Kimal Product Code	Description	French	Length (cm)	Guidewire size	Guidewire length	Arterial	Venous	Arterial	Venous	Arterial	Venous	Arterial	Medial	Venous
KFE-TDL-0608-KP	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 6Fr Length 8cm Paediatric	6	8	0.028"	50 cm	-70	90	-120	155	-190	210	0.28		0.30
KFE-TDL-0710-KP	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 7Fr Length 10cm Paediatric	7	10	0.028"	50 cm	-60	80	-100	120	-150	160	0.33		0.38
KFE-TDL-0713-KP	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 7Fr Length 13cm Paediatric	7	13	0.028"	50 cm	-70	90	-120	140	-170	190	0.39		0.41
KFE-TDL-0810-KP	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 8Fr Length 10cm Paediatric	8	10	0.028"	50 cm	-60	80	-100	120	-150	160	0.40		0.44
KFE-TDL-0813-KP	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 8Fr Length 13cm Paediatric	8	13	0.028"	50 cm	-85	100	-165	170	-220	230	0.44		0.50
KFE-TDL-0910-KP	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 9Fr Length 10cm Paediatric	9	10	0.035"	70 cm	-85	100	-165	170	-220	230	0.70		0.80
KFE-TDL-0913-KP	Short Term Dual Lumen Straight Haemodialysis Catheter Kit Size 9Fr Length 13cm Paediatric	9	13	0.035"	70 cm	-85	100	-165	170	-220	230	0.75		0.85

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

© 2024 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC.
Kimal PLC reserve the right to change the specification of its products at any time.
Kimal PLC, Arundel Road, Uxbridge, UB8 2SA.

UK
Tel: +44 (0)845 437 9542
Fax: +44 (0)845 437 9544
Email: customer care@kimal.com

International
Tel: +44 (0)845 437 9543
Fax: +44 (0)845 437 9549
Email: internationalenquiries@kimal.com